

TRACMIT User Guide

TRACMIT 1.1

<https://github.com/lacan/TRACMIT>

Contents

Installation.....	3
Dependencies	3
Using Fiji Update Sites.....	3
Manual Installation	3
Sample Dataset	3
Interface	4
ActionBar User Interface	4
Full Interface when running	5
Regions of Interest Manager	6
Parameter Optimization.....	7
Example Use	7
Bugs, problems	Error! Bookmark not defined.

Commented [P1]: Oli: please fix this.

Installation

Dependencies

As you complete the installation steps below, make sure that the following update sites are enabled:

- **IBMP-CNRS** Contains the ActionBar Plugin by Jérôme Mutterer.
- **PTBIOP** Contains the BIOPLib and attached plugins used for managing TRACMIT's settings and other internals.
- **Imagescience** Contains the FeatureJ Laplacian Plugin used by TRACMIT.

Using Fiji Update Sites

The simplest way to install TRACMIT is to use the TRACMIT Update site through Fiji:

1. From Fiji, go to **Help > Update...**
2. Select **Manage Update Sites**.
3. Click on **Add Update Site**, this will create a new line on the table.
4. Change the Name to "TRACMIT", for clarity's sake.
5. In the URL column, enter or paste <http://biop.epfl.ch/TRACMIT/>
6. Click on **Close**.
7. Finally click on **Apply Changes** and **restart Fiji**.
8. After these steps, you should find **TRACMIT** under **Plugins > ActionBar**.

Manual Installation

We do not recommend manual installation as TRACMIT depends on multiple packages that would become difficult to manage outside of the Fiji Update Site solution.

Sample Dataset

To test TRACMIT, you can download a sample dataset from ZENODO with the following DOI:

<https://doi.org/10.5281/zenodo.232218>

Interface

TRACMIT User Interface




Image Handling

Access all ASMIT Settings

Full Analysis of Current Image Stack

Run Each Step Individually
(Useful for initial parameter setting)

Inspects Currently Selected ROI

Convenience Shortcuts

Commented [P2]: This is a minor comment, which you are welcome to ignore, but perhaps it would be best to combine these two pages (4 and 5) into one, so as to have them be a bit less bulky (in comparison to the font size of the other texts)?

TRACMIT Settings Bar

 — ☐ X

All settings

Pattern

Division Detection

Mitotic Plate

Toggle Debug Mode

Wizards

FlatField Wizard

Pattern Detection Wizard

Pattern ROI Wizard


Pattern Crowding Wizard

Mitosis Detection Wizard

All settings relevant to TRACMIT.

Wizards for setting optimization

Full Interface when running




Commented [WB3]: Here it says still ASMIT. Could you change this Oli?

When TRACMIT is running, several Windows work together as shown above.

- **Original Data** contains the SIFT-Corrected original dataset loaded with Select Image or Select Raw Image.
- **Laplacian Result** is the 2D Laplacian of Gaussian timelapse used throughout TRACMIT.
- **Thresholded Laplacian Result** contains the binarized Laplacian Image used for extracting -shape features.
- **Position Mask** helps TRACMIT check if a division event was already detected at a certain location, to avoid duplicate detections.
- **Full Results** contains all the data on the detections, for each detection event at-in each frame.
- **Summary Table** is also available, with one line per division event.
- **Regions Of Interest** is the ROI Manager of ImageJ with all detections named (See Regions of Interest Manager section).
- **Settings** contains a list of all the variables used by TRACMIT. Regions of Interest Manager
- **Activity Log** helps the user follow the progress of TRACMIT.

Commented [WB4]: This should be in a separate line?

Regions of Interest Manager



At the end of processing, the ROI Manager contains 5 Types of ROIs:

- **Detected Traps:** a series of points that should match the bottom left of the patterns- ([See see](#) main publication, Figure 2 Step 5).
 - **Parameters used:** Pattern Mask Median Filter, Pattern Min Area, Pattern MaxArea.
- **Selected Candidates:** a series of points that represents the patterns kept after Standard Deviation Filtering- ([See see](#) main publication, Figure 2 Step 8).
 - **Parameters used:** Pattern Min SD, Pattern MaxSD.
- **Position Masks:** an ROI showing where the divisions will be searched for within the stack- ([See see](#) main publication, Figure 2 Step 8).
- **Potential Division #X:** Each ROI consists of two points that indicate the objects that were interpreted by TRACMIT as being anaphase figures- ([See see](#) main publication, Figure 2 Step 12). You can visit them by clicking on the **Inspect Division** Button.
 - **Parameters used:** DNA Min Area, DNA Max Area, -Division Area Tolerance, Division Angle Tolerance, Division Max Distance, ROI threshold, Max Total Area of Single Cell, Min Single Cell Area.

- **Division Event #X:** ~~Is a p~~Point ROI showing the stack position of a division event ~~in~~ which ~~where~~ a metaphase plate was found and backtracked. You can visit each division event by selecting the desired ROI and clicking on the **Inspect Division** Button- (~~See see~~ main publication, Figure 2 end of step 16)
 - **Parameters used:** Max Frame until first mitotic plate, Mitotic Plate Minimal Laplacian, Mitotic Plate Min Area, Mitotic Plate Max Area, Mitotic Plate Min Major Axis, Mitotic Plate Max Major Axis, Mitotic Plate Min Minor Axis, Mitotic Place Max Minor Axis, Mitotic Plate Min Axis Ratio, Mitotic Plate Max Axis Ratio, Mitotic Plate Max Movement, Mitotic Plate Max Frames to seek.

Parameter Optimization

Using the Parameter Wizards

The Parameter Wizards provide a simple interface to tune TRACMIT's parameters to a new dataset.

It is highly recommended to initially work on a small subset of your data. Choose regions in your data where division events are taking place and crop them out in XY and in time. This will make the initial parameter estimation and subsequent fine-tuning faster and more convenient.

In our publication, the initial datasets were 2048x2048 in size and there were 182 ~~t~~Timepoints. Cropped versions were about 512x512 ~~and across~~ 182 timepoints.

FlatField Wizard

You are prompted to draw a rectangular region ~~which that~~ corresponds to the largest structure that you would like to keep on your image. Usually, this is the size of a full micropattern if presentany, or a couple of cells long at least. This sets the **Pseudo Flatfield Blur** as the longest dimension of the drawn rectangle.

Pattern Detection Wizard

This tool ~~will~~ segments out the patterns ~~out~~ and asks s the user to select in the ROI manager (Holding the CTRL key) some representative patterns in order to estimate the pattern area ranges that TRACMIT will use to filter patterns from noise.

This sets the **Pattern Min Area** and **Pattern Max Area** -as the mean area of the selected patterns minus/plus 3 times the standard deviation, respectively.

Pattern ROI Wizard

To accommodate different patterns and geometries, this wizard averages all patterns that were detected using the criteria from the Pattern Detection Wizard.

It then prompts the user to draw a rectangle on the Average Pattern image. This pattern should represent the area inside which one expects mitosis to happen. The user is then shown what that region looks like on the detected pattern Mask.

This sets the **Division Detection Box Width**, **Division Detection Box Height**, **Division Detection Box X offset**, **Division Detection Box Y offset**, **Pattern Detection Box Width**, **Pattern Detection Box Height**

Pattern Crowding Wizard

Patterns with too many or too few cells are ~~determined-identified~~ using a Standard Deviation Projection of the time stack. With this wizard, the user is prompted to select (Holding the CTRL key) the ROIs (calculated during the Pattern Detection Wizard) that have a single cell on them. Based on this

selection, the wizard will suggest initial settings for the **Pattern Min SD** and **Pattern Max SD** as 0.7x/1.3x the computed minimum and maximum values of the selected regions, respectively.

Finally, the wizard will run the mitotic plate detection and prompt the user to instruct the wizard if the currently highlighted trap represents a crowded pattern or not. This will set the **Max Total Area of single Cell** as 1.2x the largest area the user determined as belonging to a pattern with a single cell.

Mitosis Detection Wizard

The wizard prompts the user to draw a line along the **major** axis of a metaphase plate that is going to divide in the next two frames. Using this line, the ~~Wizard~~ wizard then fits the data and calculates the metaphase plate size along the major and minor axes. It uses the fact that two frames later there is an ~~Anaphase~~ anaphase figure to estimate the acceptable Area range, the initial Laplacian of Gaussian, the Laplacian Threshold, the DNA Min and Max Areas, the Major and Minor Axis extents of the Mitotic plates and the max ~~allowed~~ movement allowed.

Settings

Laplacian Smoothing as half of the full width at half maximum of the minor axis

Laplacian Threshold as 0.2x of the minimum value of the line profile along two consecutive anaphase plates

DNA Min Area as 0.7x of the average area of each individual anaphase figure detected with the given laplacian smoothing and threshold determined above

DNA Max Area as 1.3x of the average area of each individual anaphase figure detected with the given laplacian smoothing and threshold determined above

Division Max distance as 1.2x the inter anaphase figure distance computed from the detections above

Mitotic Plate Minimal Laplacian as 1.3x the minimum value of the mitotic plate detected using the given smoothing and thresholds

Mitotic Plate Min Area as 0.7x the value of the mitotic plate area detected using the given smoothing and thresholds

Mitotic Plate Max Area as 1.3x the value of the mitotic plate area detected using the given smoothing and thresholds

Mitotic Plate Min Major Axis as 0.7x the value of the full width at half maximum of the major axis computed from the line profile.

Mitotic Plate Max Major Axis as 1.3x the value of the full width at half maximum of the major axis computed from the line profile.

Mitotic Plate Min Minor Axis as 0.7x the value of the full width at half maximum of the minor axis computed from the line profile.

Mitotic Plate Max Minor Axis as 1.3x the value of the full width at half maximum of the minor axis computed from the line profile.

Mitotic Plate Min Axis Ratio as 0.7x the ratio of the major and minor axes full width at half maximums computed from the line profile.

Mitotic Plate Max Axis Ratio as 1.3x the ratio of the major and minor axes full width at half maximums computed from the line profile.

Mitotic Plate Max Movement as 3x the xy ~~voxes-boxes~~ size

Commented [P5]: Shouldn't this be in blue and larger fonts? Also, please explain whether these are the default settings (I guess so) or recommended settings? This is not clear to me...

Commented [P6]: Correct?

Fine-Tuning the Parameters

As the interface provides a way to perform each step individually, you can optimize parameter sets independently. Refer to the section above to see which parameters are used in which step.

In order to find the values that would best match the data, such as Min/Max Areas of mitotic plates, the simplest approach is as follows:

1. Locate mitotic plates manually within the **Thresholded Laplacian Result** image.
2. Use the ImageJ Magic Wand tool to select a mitotic plate.
3. Hit “M” which will measure the current selection and give you the values associated with that mitotic plate.
4. Repeat for multiple observations and at different timepoints, and you will obtain a table ~~from which where~~ you can infer the min and max values to use.

The ~~above procedure~~ works for most parameters.

Example Use

To run our example dataset, proceed as follows:

1. Download the default parameters (TRACMIT Default Settings.txt) for this dataset from <https://github.com/lacan/TRACMIT>
2. Download the sample dataset from <https://doi.org/10.5281/zenodo.232218>.
3. Extract the dataset ZIP file.
4. Launch TRACMIT from Plugins > ActionBar > TRACMIT 1.1.
5. Click on “**Load Parameters**” and use the downloaded txt file.
6. Click on “**Select Image**”.

You will be prompted to provide the location of the extracted folder.

7. It will offer you two fields, open whichever one you would like to analyze.
8. Click on “**Measure Current Image**”.
9. The results tables might be ~~are~~ located behind an empty results table. Make sure that you move it to ~~find them see it the results table~~.
10. You can inspect the detected divisions by highlighting the one you would like to see from the ROI Manager and clicking on “**Inspect Division**”.

~~10.~~

Commented [WB7]: Could make sense to add number 11 with advice to delete any existing “processed” folders in the sample data set folder...(in case of troubles).. see e-mail.

Commented [WB8R7]:

Formatted: Font: Not Bold

Troubleshooting

Feel free to report any potential issues and to ask questions on TRACMIT's GitHub Page:

<https://github.com/lacan/TRACMIT/issues>